CLEANING

HOW TO CLEAN YOUR CAMELBAK RESERVOIR



The best way to care for your reservoir is to clean and dry it after every use, especially if you fill the reservoir with anything other than water. However, theoretically, if someone didn't clean it after *every single use*, and mold or discoloration happens to develop...

- Use hot water and 2 tablespoons of baking soda or bleach. Mix the solution inside
 your reservoir and hold it up above your head while you pinch the bite valve, allowing
 the bleached water to run through the tube. You can also use CamelBak Cleaning
 Tabs
- Let the reservoir and cleaning solution sit for about 30 minutes.
- Wash the reservoir with hot water and mild soap. Be sure to completely rinse away any bleach or cleaning solution before using again. You can also use CamelBak brushes from our Cleaning Kit to scrub your reservoir and tube. Brushes are the best way to ensure you are scrubbing all of the areas of the reservoir clean.
- Once the reservoir is clean, be sure to air dry the reservoir so no moisture is trapped inside, which can cause mold to grow.

These steps will make your reservoir safe for use. If you have spots left from the mold, they will never come all the way out because the reservoir is permanently stained, but still safe and usable.

For more information, watch our video on reservoir cleaning.

http://www.camelbak.com/Sports-Recreation/Support/Product-Care/Cleaning.aspx

Further information about Camelbak products and care can be found at:

www.camelbak.com

Jim Shaw

Lt Col, CAP

GAWG Director of Safety